

Order form LEVO C³ 2009/01

Please answer the following questions exactly and consider the indications. The exact measurement within the pictograms allows for a correct configuration and facilitates the final adjustment with the customer.
 = please mark = Standard configuration = following note



Total height: _____ cm	Diagnosis	<input type="checkbox"/> Paraplegia	<input type="checkbox"/> Quadriplegia	<input type="checkbox"/> Hemiplegia
Weight: _____ kg		<input type="checkbox"/> Multiple Sclerosis	<input type="checkbox"/> Muscl. Dystrophy	<input type="checkbox"/> Cerebral Palsy
Date of birth: _____		<input type="checkbox"/> Other: _____		
Sex: _____	F / M			

Reasons for purchasing a LEVO

Practical / Functional

- Only for standing Use as primary wheelchair
 Use at home Use as secondary wheelchair
 Use at work / school
 Other: _____

Medical / Therapeutical

- Decubitus Contractures Spasms
 Digestion Bladder problems Osteoporosis
 Circulation Pain (postural pain)
 Other: _____

How did you first hear of LEVO?

- Dealer Therapie Friends/known
 Advertisement Doctor Internet
 Other: _____

Financing


- Accident Insurance Old-age Insurance Work Insurance
 Social Security Private Insurance Charity
 Self paid Name: _____

Standard Configuration

- Measures Power Base: Width 63 cm; Length 105 cm (incl. Footplate)
- E/E Wheelchair with 4 active wheels and "Low shearing system"
- Seat support plate
- Quick release for seat depth, lower leg length and back angle
- Seat height 46 cm
- V-Trak Axxis back
- Flip-up armrests
- VR2 control system of PG-Drives Technology with 55 Ah battery
- Flip-up footrest to the side and backwards with non slip surface
- Knee support "pro"
- Vehicle tie-down locations

1. Upper leg length

S - XL



Measure while sitting the distance between the rear point of the bottom and the hollow of the knee.

38 cm ← → 66 cm



Adjustable with quick release over the hole range.



Attention: Do not confound with the seat depth which is 3 cm shorter!

2. Hip size



Measure in sitting position the outermost points from the bottom or thigh.

32 cm ← → 52 cm

Seat width:

- | | | |
|-------|---------------|--------------------------|
| 32 cm | (Art. Number) | <input type="checkbox"/> |
| 36 cm | (Art. Number) | <input type="checkbox"/> |
| 40 cm | (Art. Number) | <input type="checkbox"/> |
| 44 cm | (Art. Number) | <input type="checkbox"/> |
| 48 cm | (Art. Number) | <input type="checkbox"/> |
| 52 cm | (Art. Number) | <input type="checkbox"/> |

with chest strap "Bodypoint" 150 cm long / 6 cm wide

(Art. Number)

(Standard for seat width 32-40 cm)

with chest strap "Bodypoint" 180 cm long / 8 cm wide

(Art. Number)

(Standard for seat width 44-52 cm)

3. Lower leg length



Please consider the height of the cushion.
 The LEVO cushion increases the distance from the seat to the footplate by 6 cm. For measuring **Shoes** should be worn.

24 cm ← → 44 cm

Option:

With lower leg length between 45 and 50 cm the wheelchair is available with 6 cm higher total seat height. The seat height is fix and not adjustable.

(Art. Number)



Calculation example:

lower leg length	49 cm
minus cushion height	-6 cm
Footrest position	43 cm

4. Backrest



V-Trak Axxis back is used as a standard. Of course other back systems can be used.

V-Trak Axxis Reduced Height Backrest, 31 cm/12" high
 V-Trak Axxis Low Backrest, 40 cm/16" high
 V-Trak Axxis High Backrest, 61 cm/24" high

(Art. Number)	<input type="checkbox"/>
(Art. Number)	<input type="checkbox"/>
(Upcharge)	(Art. Number) <input type="checkbox"/>

without back system

5. Armrest height (measured with LEVO cushion 6cm high)

15 cm ← → 26 cm ◆ 22 cm

6. Operation (Joystick modul)

◆ right left

7. Max. Speed

6 km/h 8 km/h ◆ 10 km/h

8. Driving characteristics

Classic ◆ Standard Fun

9. Battery capacity

range approx. 25 km 55 Ah ◆ (Art. Number)
 range approx. 35 km 72 Ah (Upcharge) (Art. Number)

10. Frame colors

Other: RAL-Nr.: _____ (Upcharge) (Art. Number)

red (Art. Number)
 yellow (Art. Number)
 blue (Art. Number)
 black (Art. Number)

red metallise (Upcharge) (Art. Number)
 green metallise (Upcharge) (Art. Number)
 blue metallise (Upcharge) (Art. Number)
 silver metallise (Upcharge) (Art. Number)

orange reflex (Upcharge) (Art. Number)
 lemon reflex (Upcharge) (Art. Number)
 silver reflex (Upcharge) (Art. Number)

11. Options

Knee support



◆ Knee support "pro"

Foam size S (10-12 cm)
 Foam size M (14-16 cm)

Knee support "integral"



(Upcharge) (Art. Number)

Foam size S (10-12 cm)
 Foam size M (14-16 cm)

Swing away joystick



- Creates space for transfer
 - Possible to drive closer to tables

(Upcharge) (Art. Number)

Joystick protection guard



Shock protection for the joystick. Protects also the hands of the user. **Only possible when ordering with VR2!**

(Upcharge) (Art. Number)



R-Net Control



- 120 A
 - High resolution monochrome graphics
 LCD Display Connector block: For additional Devices (Dual Attendant, environmental control)

(Price) (Art. Number)



Colour display



Using R-net, connector points for accessory "Buddy Button" are included. With accessory "Omni Module" it is possible to connect environmental controls of different supplier.

(Upcharge) (Art. Number)



Omni module



- Compatible with many input devices
- Remote control
- 4/5 switch control
- 3 axis switch or proportional (head control)
- Sip and puff
- Single switch scanner

(Price) (Art. Number)



Arms to adapt the Omni



Flexible arm incl. Adapter

(Price) (Art. Number)

Joint arm incl. Adapter

(Price) (Art. Number)



Chin control with shoulder holder



Shoulder- and chest holder with upholstery. Height and width individually adjustable. With "ON/OFF" und "MODE" button.

(Price) (Art. Number)



Chin control to deviate SPZ



Fixed at the chair and electrical deviation. Adjustable in height, width and angle.

(Price) (Art. Number)



Dual attendant module




Actuation for assisting people. Easy handling and switch over.

(Price) (Art. Number)



Lamps for outdoor use with indicators



(Upcharge) (Art. Number)

Tilt in Space (0 - 35°)



- comfortable relax position
- prevention pressure sores

(Upcharge) (Art. Number)

Footplate, depth 20 cm



Standard configuration is 25 cm deep. This option reduces the total chair length.

(Upcharge) (Art. Number)

Wheels



V-Profile for drive wheels



For more indoor use.

(Art. Number)



Block-profile for drive wheels



For more outdoor use.

(Art. Number)

12. Accessories

Foam seat cushion with seat cushion cover



Seat cushion corresponding to upper leg length and hip size

(Price) (Art. Number)

additional seat cushion cover

(Price) (Art. Number)

Stroller handle



Can be used with standard and V-Trak back. Helps for manual pushing of the wheelchair, e.g. to load into a vehicle.

(Price) (Art. Number)

Swing away table tray



(Price) (Art. Number)

Skirt guard



Easy to remove for lateral transfer.

Pair (Price) (Art. Nombre)

Upper leg support



Can be adapted individually in place, distance, angle and height. Leads the thigh and stabilizes the leg axle in seating and standing

Pair (Price) (Art. Number)

Chest role



Individually adjustable in height. Gives a feeling of more safety and stability, especially while standing. Does not replace the chest strap.

(Price) (Art. Number)

Foot guide



Can be adjusted individually. Allows a better foot position and stability.

small right
small left

(Price) (Art. Number)
(Price) (Art. Number)

large right
large left

(Price) (Art. Number)
(Price) (Art. Number)

Buddy Button



green: ON/OFF switch for joystick module
yellow: MODE switch for joystick module
red: Emergency stop, disconnect the power supply

green
yellow
red

(Price) (Art. Number)
(Price) (Art. Number)
(Price) (Art. Number)

Joystick-top



T-bar
Stick
Softball
Mushroom
Chin
Ball

(Price) (Art. Number)
(Price) (Art. Number)
(Price) (Art. Number)
(Price) (Art. Number)
(Price) (Art. Number)
(Price) (Art. Number)

Rear view mirror



Joystick side
opposite joystick

(Price)	(Art. Number)	<input type="checkbox"/>
(Price)	(Art. Number)	<input type="checkbox"/>

Monitor



For an objective outcome measurement for users and therapists. Delivered with PC-Software.

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

At option V-Track Headrest (incl. Track extension):



V-Trak Headrest winged
V-Trak Headrest simple

(Price)	(Art. Number)	<input type="checkbox"/>
(Price)	(Art. Number)	<input type="checkbox"/>

Additional information (ev. diagram):

Total: _____

Order No:

Date:

LEVO-Dealer:

Client: