

# Order form **LEVO LAE** 2009/01



Please answer the following questions exactly and consider the indications. The exact measurement within the pictograms allows for a correct configuration and facilitates the final adjustment with the customer.

= please mark  = Standard configuration ⓘ = following note

Total height: \_\_\_\_\_ cm     **Diagnosis**      Paraplegia      Quadriplegia      Hemiplegia  
 Weight: \_\_\_\_\_ kg      Multiple Sclerosis      Muscl. Dystrophy      Cerebral Palsy  
 Date of birth: \_\_\_\_\_      Other: \_\_\_\_\_  
 Sex:     F / M     \_\_\_\_\_

## Reasons for purchasing a LEVO


<b>Practical / Functional</b> <input type="checkbox"/> Only for standing <input type="checkbox"/> Use as primary wheelchair <input type="checkbox"/> Use at home <input type="checkbox"/> Use as secondary wheelchair <input type="checkbox"/> Use at work / school <input type="checkbox"/> Other: _____	<b>Medical / Therapeutical</b> <input type="checkbox"/> Decubitus <input type="checkbox"/> Contractures <input type="checkbox"/> Spasms <input type="checkbox"/> Digestion <input type="checkbox"/> Bladder problems <input type="checkbox"/> Osteoporosis <input type="checkbox"/> Circulation <input type="checkbox"/> Pain (postural pain) <input type="checkbox"/> Other: _____
---	---

<b>How did you first hear of LEVO?</b> <input type="checkbox"/> Dealer <input type="checkbox"/> Therapie <input type="checkbox"/> Friends/known <input type="checkbox"/> Advertisement <input type="checkbox"/> Doctor <input type="checkbox"/> Internet <input type="checkbox"/> Other: _____	<b>Financing</b> <input type="checkbox"/> Accident Insurance <input type="checkbox"/> Old-age Insurance <input type="checkbox"/> Work Insurance <input type="checkbox"/> Social Security <input type="checkbox"/> Private Insurance <input type="checkbox"/> Charity <input type="checkbox"/> Pay self <input type="checkbox"/> Name: _____
---	--

## Standard Configuration

- Seat- and back upholstery
- One arm drive right or left
- Two piece flip up footrest
- Pneumatic rear wheels
- Fold down backrest / angle adjustable
- Chest strap and knee support

### 1. Upper leg length



Measure in sitting position the distance between the rear point of the bottom and the hollow of the knee.

**S - M** (adjustable for upper leg length 35 - 48 cm)

35 cm ← → 48 cm

◆ S = 35 cm  
 ◆ M = 48 cm

(Art. Number)   
 (Art. Number)

**L - XL** (adjustable for upper leg length 48 - 63 cm)

49 cm ← → 63 cm

◆ L = 52 cm  
 ◆ XL = 63 cm

(Art. Number)   
 (Art. Number)



**Attention:** Do not confound with the seat depth which is 6-8 cm shorter !

### 2. Hip size



Measure in sitting position the outermost points from the bottom or thigh.

35 cm ← → 48 cm

#### Seat width:

35 cm	(Art. Number) <input type="checkbox"/>
38 cm	(Art. Number) <input type="checkbox"/>
40 cm	(Art. Number) <input type="checkbox"/>
42 cm	(Art. Number) <input type="checkbox"/>
45 cm	(Art. Number) <input type="checkbox"/>
48 cm	(Upcharge) (Art. Number) <input type="checkbox"/>


with chest strap "Bodypoint" 150 cm long / 6 cm wide     (Art. Number)

(Standard for seat width 35-42 cm)

with chest strap "Bodypoint" 180 cm long / 8 cm wide     (Art. Number)

(Standard for seat width 45-48 cm)

### 3. Lower leg length



#### Please consider the height of the cushion.

The LEVO cushion increases the distance from the seat to the footplate by 3 cm / 1.2". For the measuring **Shoes** should be worn.

36 cm ← → 56 cm

Depending on the lower leg length, you can choose between two or three different seat heights. The selected seat height affects the wheel size and thus also the possibility to drive under the table, but also the user's expenditure of energy.

#### Footrest position:

Seat height 48	36 cm	(Art. Number) <input type="checkbox"/>
Seat height 51	38 cm	(Art. Number) <input type="checkbox"/>
Seat height 54	40 cm	(Art. Number) <input type="checkbox"/>
Seat height 57	42 cm	(Art. Number) <input type="checkbox"/>
	44 cm	(Art. Number) <input type="checkbox"/>
	46 cm	(Art. Number) <input type="checkbox"/>
	48 cm	(Art. Number) <input type="checkbox"/>
	50 cm	(Art. Number) <input type="checkbox"/>
	52 cm	(Art. Number) <input type="checkbox"/>
	54 cm	(Art. Number) <input type="checkbox"/>
	56 cm	(Art. Number) <input type="checkbox"/>




**Calculation example:** lower leg length 51.0 cm  
 minus cushion height -3.0 cm  
**Footrest position 48.0 cm** Choose between seat height 51/54 and 57 cm

### 4. Seat - floor height front (rear) without cushion

Rear wheel 22", castor 5"	48 (45 cm)	(Art. Number) <input type="checkbox"/>
Rear wheel 24", castor 5"	51 (47 cm)	(Art. Number) <input type="checkbox"/>
Rear wheel 24", castor 5"	54 (49 cm)	(Art. Number) <input type="checkbox"/>
Rear wheel 26", castor 6"	57 (52 cm)	(Upcharge) (Art. Number) <input type="checkbox"/>

### 5. Backrest height



Distance between seating surface and the lower blade bone.

**Standard Backrest:**

- 31 cm   
  33 cm   
  35 cm   
  37 cm   
  39 cm   
  41 cm

**Option:**

- Backrest extension 15 cm  
 V-Trak Standard Backrest, 40 cm/16" high  
 V-Trak Medium Backrest, 60 cm/24" high \*  
 V-Trak High Backrest, 80 cm/31.5" high \*  
 V-Trak Axxis Reduced Height Backrest, 31 cm/12" high  
 V-Trak Axxis Low Backrest, 40 cm/16" high  
 V-Trak Axxis High Backrest, 61 cm/24" high  
 V-Trak Backrest extension, 20 cm  
 \* Track extension at headrest use

(Price)	(Art. Number)	<input type="checkbox"/>
(Upcharge)	(Art. Number)	<input type="checkbox"/>
(Upcharge)	(Art. Number)	<input type="checkbox"/>
(Upcharge)	(Art. Number)	<input type="checkbox"/>
(Upcharge)	(Art. Number)	<input type="checkbox"/>
(Upcharge)	(Art. Number)	<input type="checkbox"/>
(Upcharge)	(Art. Number)	<input type="checkbox"/>
(Price)	(Art. Number)	<input type="checkbox"/>

### 6. Backrest angle

- 85°                      
  90°                      
 95°

### 7. Stand-up operation




Since 2007 only one-hand operation available!

- right                     
  left

### 8. Footplate angle (as a standard adjustable into all three positions)

- 80°                      
  90°                      
 100°

### 9. Rear wheel position



- Position S                      
 Position M                      
 Position L                      
 Position XL



The shorter the wheel base is, the more dynamical is the handling of the wheel-chair, but also the tilt risk is increased. For example, you can choose for wheelchair size "L" the rear wheel position "M". In this case anti-tip aid are recommended.

### 10. Gas springs



The stronger the gas springs are, the more support you have when standing up, but also more resistance to sit down.

- For body weight: 40-60 kg/88-132 pounds                      400 N  (Art. Number)  
 For body weight: 60-80 kg/132-176 pounds                      500 N  (Art. Number)  
 (Standard for seat depth M)  
 For body weight: 80-100 kg/176-220 pounds                      600 N  (Art. Number)  
 (Standard for seat depth L)  
 For body weight: 100-120 kg/220-264 pounds                      700 N  (Art. Number)  
 (Standard for seat depth XL)

### 11. Frame colors

- platinum  (Art. Number)                     
 Other: RAL-No.: \_\_\_\_\_  (Upcharge)  (Art. Number)  
 red  (Art. Number)                     
 red metallise  (Upcharge)  (Art. Number)  
 yellow  (Art. Number)                     
 green metallise  (Upcharge)  (Art. Number)                     
 orange reflex  (Upcharge)  (Art. Number)  
 blue  (Art. Number)                     
 blue metallise  (Upcharge)  (Art. Number)                     
 lemon reflex  (Upcharge)  (Art. Number)  
 black  (Art. Number)                     
 silver metallise  (Upcharge)  (Art. Number)                     
 silver reflex  (Upcharge)  (Art. Number)  
 aluminium polished  (Upcharge)  (Art. Number)

### 12. Knee wedge



Knee supports are only delivered extra strong. Stabilization is removable.

#### Rectangular (Standard)



- 50x100  (Art. Number)  
 60x100  (Art. Number)  
 80x100  (Art. Number)  
 50x150  (Art. Number)  
 80x150  (Art. Number)

#### Conical (Option)



The conical knee wedges give more support, especially with spasticity.

- 50-100x100  (Upcharge)  (Art. Number)  
 90-140x100  (Upcharge)  (Art. Number)  
 130-180x100  (Upcharge)  (Art. Number)

### 13. Options

#### Armrest colors




















- black  (Art. Number)                     
 Other: RAL-No.: \_\_\_\_\_  (Upcharge)  (Art. Number)  
 platinum  (Upcharge)  (Art. Number)                     
 red metallise  (Upcharge)  (Art. Number)                     
 orange reflex  (Upcharge)  (Art. Number)  
 red  (Upcharge)  (Art. Number)                     
 green metallise  (Upcharge)  (Art. Number)                     
 lemon reflex  (Upcharge)  (Art. Number)  
 yellow  (Upcharge)  (Art. Number)                     
 blue metallise  (Upcharge)  (Art. Number)                     
 silver reflex  (Upcharge)  (Art. Number)  
 blue  (Upcharge)  (Art. Number)                     
 silver metallise  (Upcharge)  (Art. Number)                     
 aluminium polished  (Upcharge)  (Art. Number)

#### Legrest removable for transport



The benefit of being able to remove leg rest by quick release is, to reduce the volume for transportation.

- (Upcharge)  (Art. Number)

<b>One-piece footrest fix, no flip-up feature</b> 			<table border="1"> <tr> <td>-----</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> </table>	-----	(Art. Number)	<input type="checkbox"/>						
-----	(Art. Number)	<input type="checkbox"/>										
<b>Ground level footrest</b> 		Especially for users who use the wheelchair as a stand-up help.	<table border="1"> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> </table>	(Upcharge)	(Art. Number)	<input type="checkbox"/>						
(Upcharge)	(Art. Number)	<input type="checkbox"/>										
<b>Seat depth / Footrest quick adjustment</b> (Set: 10 Clips, 4 Screws) 		Is suitable for demo chairs especially, which are used for different clients.	<table border="1"> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> </table>	(Upcharge)	(Art. Number)	<input type="checkbox"/>						
(Upcharge)	(Art. Number)	<input type="checkbox"/>										
<b>Footrest height adapter</b> 		SB: 35, 38, 45, 48 cm SB: 40, 42 cm *  Especially for users with short lower leg length. The footrest becomes 3 cm narrower than the seat width.  * <b>Only</b> possible with "one-piece footplate" .	<table border="1"> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> </table>	(Upcharge)	(Art. Number)	<input type="checkbox"/>	(Upcharge)	(Art. Number)	<input type="checkbox"/>			
(Upcharge)	(Art. Number)	<input type="checkbox"/>										
(Upcharge)	(Art. Number)	<input type="checkbox"/>										
<b>Non-Slip surface</b>		Available for standard footrest (two piece flip-up), but also for one piece footrest fix.	<table border="1"> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> </table>	(Upcharge)	(Art. Number)	<input type="checkbox"/>						
(Upcharge)	(Art. Number)	<input type="checkbox"/>										
<b>Rear wheels with airless tyres</b>		Standard: Air tyres with aluminium hand ring.	<table border="1"> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> </table>	(Upcharge)	(Art. Number)	<input type="checkbox"/>						
(Upcharge)	(Art. Number)	<input type="checkbox"/>										
<b>Rear wheels with hub locking brake</b> 		Operation without additional energy. It does not depend on the tyre air pressure. Less frame weight for transportation.	<table border="1"> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> </table>	(Upcharge)	(Art. Number)	<input type="checkbox"/>						
(Upcharge)	(Art. Number)	<input type="checkbox"/>										
<b>Rear wheels with drum breaks</b> 		Possibility for assistants to stop the wheelchair, safe and comfortable.	<table border="1"> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> </table>	(Upcharge)	(Art. Number)	<input type="checkbox"/>						
(Upcharge)	(Art. Number)	<input type="checkbox"/>										
<b>Rear wheels "Spinergy"</b> 		Light weight, improved driving performance, easy to clean and better design.	<table border="1"> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> </table>	(Upcharge)	(Art. Number)	<input type="checkbox"/>	(Upcharge)	(Art. Number)	<input type="checkbox"/>	(Upcharge)	(Art. Number)	<input type="checkbox"/>
(Upcharge)	(Art. Number)	<input type="checkbox"/>										
(Upcharge)	(Art. Number)	<input type="checkbox"/>										
(Upcharge)	(Art. Number)	<input type="checkbox"/>										
<b>Camber</b> 		Improved driving performance and agility. Totally width is slightly increased.	<table border="1"> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> </table>	(Upcharge)	(Art. Number)	<input type="checkbox"/>	(Upcharge)	(Art. Number)	<input type="checkbox"/>			
(Upcharge)	(Art. Number)	<input type="checkbox"/>										
(Upcharge)	(Art. Number)	<input type="checkbox"/>										
<b>Steel hand ring</b> 		Increased grip, less temperature-sensitive and less dirty hands.	<table border="1"> <tr> <td>(Upcharge)</td> <td>(Art. Number)</td> <td><input type="checkbox"/></td> </tr> </table>	(Upcharge)	(Art. Number)	<input type="checkbox"/>						
(Upcharge)	(Art. Number)	<input type="checkbox"/>										

**Maxgrepp** (Hand ring with rubber covered)



**Original**

**Ergo Para CP**



Better hold for persons with limited strength or grip.

Pair

(Upcharge)	(Art. Number)	<input type="checkbox"/>
------------	---------------	--------------------------

Pair

(Upcharge)	(Art. Number)	<input type="checkbox"/>
------------	---------------	--------------------------

**One arm drive**



For stroke patients or people with unilateral

right  
left

(Upcharge)	(Art. Number)	<input type="checkbox"/>
------------	---------------	--------------------------

(Upcharge)	(Art. Number)	<input type="checkbox"/>
------------	---------------	--------------------------



Only possible with standard wheel.

**14. Accessory**

**Foam seat cushion with seat cushion cover**



S 29 cm depth  
M 38 cm depth  
L 43 cm depth  
XL 52 cm depth

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

additional seat cushion cover

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

**Push handle**



**Push handle**

Pair

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------



**Push handle extension adjustable**

Pair

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------



Helpful if the wheelchair is pushed by different persons.

**Anti-tip**



**Anti-tip rolls**

right  
left

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------



**Anti-tip aid / tip aid swing away "HC"**

right  
left

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------



**XL adapter**  
(needed for rear wheel position XL)

right  
left

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------



**Tip aid**

right  
left

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------

**Foot guide**



Individually adjustable. Allows a better foot position and foot stability.

small right  
small left

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

large right  
large left

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

**Brake lever extension (6 cm/2.4")**



Simpler accessibility and less strength necessary.

right  
left

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

**Spoke guard**



**Spoke guard "Transparent"**

right  
left

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

**Spoke guard "Transparent with LEVO-Logo"**

right  
left

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

**Vehicle tie-down kit**



Fixed attachment points for chair transport in the car. Must be combined with frame reinforcement.

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

**Frame reinforcement for e-fix or the transportation of persons in wheelchair**



Wheelchairs with e-fix need frame reinforcement. You can only have this combined with anti-tip aid "HC" and **not** with anti-tip rolls or tip aid.

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

**Monitor**



For an objective outcome measurement for users and therapists. Delivered with PC-Software.

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

**At option V-Trak:**

**V-Trak Thoracic support & pad**



- V-Trak Thoracic support & pad & cover large right
- V-Trak Thoracic support & pad & cover large left
- V-Trak Thoracic support & pad & cover medium right
- V-Trak Thoracic support & pad & cover medium left

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------

**V-Trak Headrest (incl. Trak extension)**



- V-Trak Headrest winged
- V-Trak Headrest simple

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

**Additional information** (ev. diagram):

**Total:** \_\_\_\_\_

**Order No:** \_\_\_\_\_

**Date:** \_\_\_\_\_

**LEVO-Dealer:** \_\_\_\_\_

**Client:** \_\_\_\_\_