

Order form **LEVO** kid&junior 2009/01



Please answer the following questions exactly and consider the indications. The exact measurement within the pictograms allows for a correct configuration and facilitates the final adjustment with the customer.

= please mark = Standard configuration = following note

Total height: _____ cm **Diagnosis** Paraplegia Quadriplegia Hemiplegia
 Weight: _____ kg Multiple Sclerosis Muscl. Dystrophy Cerebral Palsy
 Date of birth: _____
 Sex: F / M

Reasons for purchasing a LEVO


Practical / Functional <input type="checkbox"/> Only for standing <input type="checkbox"/> Use as primary wheelchair <input type="checkbox"/> Use at home <input type="checkbox"/> Use as secondary wheelchair <input type="checkbox"/> Use at work / school <input type="checkbox"/> Other: _____	Medical / Therapeutical <input type="checkbox"/> Decubitus <input type="checkbox"/> Contractures <input type="checkbox"/> Spasms <input type="checkbox"/> Digestion <input type="checkbox"/> Bladder problems <input type="checkbox"/> Osteoporosis <input type="checkbox"/> Circulation <input type="checkbox"/> Pain (postural pain) <input type="checkbox"/> Other: _____
---	---

How did you first hear of LEVO? <input type="checkbox"/> Dealer <input type="checkbox"/> Therapie <input type="checkbox"/> Friends/Known <input type="checkbox"/> Advertisement <input type="checkbox"/> Doctor <input type="checkbox"/> Internet <input type="checkbox"/> Other: _____	Financing <input type="checkbox"/> Accident Insurance <input type="checkbox"/> Old-age Insurance <input type="checkbox"/> Work Insurance <input type="checkbox"/> Social Security <input type="checkbox"/> Private Insurance <input type="checkbox"/> Charity <input type="checkbox"/> Pay self <input type="checkbox"/> Name: _____
---	--

Standard Configuration

- Seat- and back upholstery
- Armrests folding upwards
- Armrests height and angle adjustable
- Two piece flip up footrest
- Pneumatic rear wheels
- Fold down backrest / angle adjustable
- Battery and charger
- Chest strap and knee support

1. Upper leg length



Measure in sitting position the distance between the rear point of the bottom and the hollow of the knee.

XS - S (adjustable for upper leg length 28 - 45 cm)

28 cm ← _____ → 45 cm


XS = 24 cm
 S = 48 cm

(Art. Number)
(Art. Number)



Attention: Do not confound with the seat depth which is 6-8 cm shorter !

2. Hip size




Measure in sitting position the outermost points from the bottom or thigh.

< 30 cm ← _____ → 35 cm

Seat width:

30 cm	(Art. Number)	<input type="checkbox"/>
35 cm	(Art. Number)	<input type="checkbox"/>

3. Lower leg length



Please consider the height of the cushion. The Levo cushion increases the distance from the seat to the footplate by 2.5 cm/1". For the measuring **Shoes** should be worn.

19 cm ← _____ → 41 cm

Footrest position:

Kid 33 with Adapter	19.2 cm	(Art. Number)	<input type="checkbox"/>
Jun. w. o. Ad.	21.0 cm	(Art. Number)	<input type="checkbox"/>
Junior 40 with Adapter	22.8 cm	(Art. Number)	<input type="checkbox"/>
	24.6 cm	(Art. Number)	<input type="checkbox"/>
	26.4 cm	(Art. Number)	<input type="checkbox"/>
Junior 44 w. o. Adapter	28.2 cm	(Art. Number)	<input type="checkbox"/>
Junior 44 with Adapter	30.0 cm	(Art. Number)	<input type="checkbox"/>
	31.8 cm	(Art. Number)	<input type="checkbox"/>
	33.6 cm	(Art. Number)	<input type="checkbox"/>
	35.4 cm	(Art. Number)	<input type="checkbox"/>
	37.2 cm	(Art. Number)	<input type="checkbox"/>
	39.0 cm	(Art. Number)	<input type="checkbox"/>
	40.8 cm	(Art. Number)	<input type="checkbox"/>

Footrest adapter for short legs (Standard with *kid*)

(Up charge)	(Art. Number)	<input type="checkbox"/>
-------------	---------------	--------------------------




Calculation example: lower leg length 36.0 cm
 minus cushion height -2.5 cm
Footrest position 33.5 cm Choose between seat height 33/40 or 44 cm

4. Seat - floor height front (rear) without cushion

Rear wheel 16", castor 3"	<i>kid</i> : 33 (30) cm	(Art. Number)	<input type="checkbox"/>
Rear wheel 20", castor 5"	<i>junior</i> : 40 (37) cm	(Art. Number)	<input type="checkbox"/>
Rear wheel 22", castor 5"	<i>junior</i> : 44 (41) cm	(Art. Number)	<input type="checkbox"/>

5. Backrest height



Distance between seating surface and the lower blade bone.

Standard Backrest:

◆ 28 cm

(Art. Number)

with backrest extension:

29 cm



42 cm

(Up charge)

(Art. Number)

Option:

V-Trak 28 cm/11" (only for seat width 30 cm)

(Up charge) (Art. Number)

V-Trak 30 cm/12"

(Up charge) (Art. Number)

V-Trak Axxis Reduced Backrest Height 31 cm/12"

(Up charge) (Art. Number)

V-Trak Axxis Low Backrest 40 cm/16"

(Up charge) (Art. Number)

V-Trak Backrest extension +13 cm

(Up charge) (Art. Number)

6. Backrest angle

85°

◆ 90°

95°

7. Stand-up mechanism



Standard on standing operation (armrest round with control switch).

◆ right

left

8. Armrest pad (measured without cushion)



10 cm



22 cm

(stepless height adjustable)

9. Footplate angle

70°

80°

◆ 90°

100°

110°

10. Frame colors

Andere: RAL-Nr.:

(Up charge) (Art. Number)

platin (Art. Number)
 red (Art. Number)
 yellow (Art. Number)
 blue (Art. Number)
 black (Art. Number)

red métallisé (Up charge) (Art. Number)
 green métallisé (Up charge) (Art. Number)
 blue métallisé (Up charge) (Art. Number)
 silver métallisé (Up charge) (Art. Number)

orange reflex (Up charge) (Art. Number)
 lemon reflex (Up charge) (Art. Number)
 silverreflex (Up charge) (Art. Number)
 aluminium polished (Up charge) (Art. Number)

11. Knee wedge



Knee supports are only delivered extra strong. Stabilization is removable.

Rectangular (Standard)



50x100 (Art. Number)
 60x100 (Art. Number)
 50x150 (Art. Number)

Conical (Option)



The conical knee wedges give more support, especially with spasticity.

50-100x100 (Up charge) (Art. Number)
 90-140x100 (Up charge) (Art. Number)
 130-180x100 (Up charge) (Art. Number)

12. Options

Dual switch for table control



Stand-up operation with push buttons for table. Table tray not included (see accessories).

(Up charge) (Art. Number)

Hand control



Instead of control switch

-. (Art. Number)

Seat depth / Footrest quick adjustment (Set: 10 Clips with Bolt, 4 Screws)



Is suitable for demo chairs especially, which are used for different clients.

(Up charge) (Art. Number)

Non-Slip surface

(Up charge) (Art. Number)

Rear wheels "Spinergy"



Only possible with 22"-wheels.

red
black
white

(Up charge) (Art. Number)
 (Up charge) (Art. Number)
 (Up charge) (Art. Number)



Light weight, improved driving performance, easy to clean and better design.

Maxgrepp (Hand ring with rubber covered)



Original



Only possible with 22"-wheels.

Pair

(Up charge) (Art. Number)



Ergo Para CP



Better hold for persons with limited strenght or

Pair

(Up charge) (Art. Number)

13. Accessory

Foam seat cushion with seat cushion cover



XS 30	length 25 x width 30	(Price)	(Art. Number)	<input type="checkbox"/>
XS 35	length 25 x width 35	(Price)	(Art. Number)	<input type="checkbox"/>
S 30	length 34 x width 30	(Price)	(Art. Number)	<input type="checkbox"/>
S 35	length 35 x width 35	(Price)	(Art. Number)	<input type="checkbox"/>
additional seat cushion cover		(Price)	(Art. Number)	<input type="checkbox"/>

Table tray



(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

Control switch extension



vertical	(Price)	(Art. Number)	<input type="checkbox"/>
horizontal	(Price)	(Art. Number)	<input type="checkbox"/>

Height adjustable stroller handle



30 cm	(Price)	(Art. Number)	<input type="checkbox"/>	
35 cm	(Price)	(Art. Number)	<input type="checkbox"/>	
with stand-up operation for attendant	30 cm	(Price)	(Art. Number)	<input type="checkbox"/>
	35 cm	(Price)	(Art. Number)	<input type="checkbox"/>

Folding down stroller bar HC (width adjustable)



(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

Anti-tip roll



Anti-tip for kid	(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
-------------------------	---------	------------	---------------	--------------------------

Anti-tip for jun	(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
-------------------------	---------	------------	---------------	--------------------------

Back pack



red	(Price)	(Art. Number)	<input type="checkbox"/>
yellow	(Price)	(Art. Number)	<input type="checkbox"/>
dark blue	(Price)	(Art. Number)	<input type="checkbox"/>
black	(Price)	(Art. Number)	<input type="checkbox"/>

Foot guide



Individually adjustable. Allows a better foot position and foot stability.

small right	(Price)	(Art. Number)	<input type="checkbox"/>
small left	(Price)	(Art. Number)	<input type="checkbox"/>

Knee support height adapter



Is suitable for persons with long lower legs. Knee holder gets closer to the joint, better pressure distribution.

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

Brake lever extension (6 cm/2.4")



Easier accessibility and less strength necessary.

right	(Price)	(Art. Number)	<input type="checkbox"/>
left	(Price)	(Art. Number)	<input type="checkbox"/>

Spoke guard



Spoke guard "Tiger"



Only possible with 20" and 22"-wheels

Pair	(Price)	(Art. Number)	<input type="checkbox"/>
------	---------	---------------	--------------------------

Monitor



For an objective outcome measurement for users and therapists. Delivered with PC-Soft-ware.

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

At option V-Trak:

V-Trak Thoracic support & pad



- V-Trak Thoracic support & pad & cover large right
- V-Trak Thoracic support & pad & cover large left
- V-Trak Thoracic support & pad & cover medium right
- V-Trak Thoracic support & pad & cover medium left

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------

(Price)	(Quantity)	(Art. Number)	<input type="checkbox"/>
---------	------------	---------------	--------------------------

V-Trak Headrest (incl. Trak extension)



- V-Trak Headrest winged
- V-Trak Headrest simple

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

(Price)	(Art. Number)	<input type="checkbox"/>
---------	---------------	--------------------------

Additional information (ev. diagram):

Total: _____

Order No:

Date:

LEVO-Dealer:

Client: