

LEVO Summit EL

The next generation of the LEVO manual drive standing wheelchairs! Summit EL – powered standing as exceptional performance and style!

LEVO Summit EL



LEVO stands for a better quality of life

LEVO[®].dt 
THE EXPERTS IN STANDING

Inspired Swiss design

Innovative “ABC” standing

The first truly active standing wheelchair

This is LEVO: combining 50 years of development and production experiences, LEVO is proud to offer “best in class” manual standing wheelchair line.





Our innovative design delivers the features of an active wheelchair and state of the art standing technology, allowing you to reach the perfect standing position, from any seating position you choose.

The **LEVO Summit** combines the lowest stand-up wheelchair weight with the highest user capacity in the market.

Highlights:

- Weight from 21 kg
- User weight capacity max. 125 kg
- Automatic backrest compensation “ABC” for ergo correct standing
- Easily configurable and convertible to individual needs
- Converts from manual/manual to manual/electric
- Easy transfer



Technical Data

Model	S - L
Seat width (SW)	35 / 38 / 40 / 42 / 45 / 48 / 50
Total width (W)	Seat width + 16 cm
Total length (L)	Seat depth + 43 cm
Total height (H) min.	Seat height + 30 cm
Backrest upholstery height	32 – 38 cm
Seat height	48 / 50 / 52 / 54
Seat depth	35 – 55 with 2 cm increments
Footrest	height and angle adjustable
Total weight	min. 21 kg
Batteries	24 Volt, 2.9 Ah Li-Ion

Distributor:

Manufactured by:

LEVO AG | Anglikerstrasse 20 | Tel +41 (0)56 618 44 11 | office@levo.ch | UID-Nr. CHE-107.867.006 MWST
Switzerland | CH-5610 Wohlen | Fax +41 (0)56 618 44 10 | www.levo.ch | EU SRN CH-MF-000008751

LEVO.^{dt}
THE EXPERTS IN STANDING