

LEVO Summit EL

The next generation of the LEVO LCEV – Delivering exceptional performance and style

LEVO Summit EL



LEVO stands for a better quality of life

LEVO.dt⁺
THE EXPERTS IN STANDING

Inspired Swiss design

Innovative “ABC” standing

The first truly active standing wheelchair

This is LEVO: combining 40 years of market leading standing experience LEVO is proud to introduce its new “best in class” manual standing wheelchair line.




Our new innovative design delivers the features of an active wheelchair and state of the art standing technology, allowing you to reach the perfect standing position, from any seating position you choose.

The **LEVO Summit** combines the lowest stand-up wheelchair weight with the highest user capacity in the market.

Highlights:

- Weight of just 46.3 lbs
- Weight capacity max. 275 lbs
- Automatic backrest compensation “ABC” for ergo correct standing
- Compact active driving
- Easily configurable and convertible to individual needs
- Converts from manual/manual to manual/electric
- Easy transfer
- Easy transport



Technical Data

Model	S-L
Seat width (SW)	13.8 / 15 / 15.7 / 16.5 / 17.7 / 18.9 / 19.7
Total width (W)	Seat width + 6.3 in
Total length (L)	Seat depth + 16.9 in
Total height (H) min.	Seat high + 11.8 in
Backrest upholstery height	12.6 – 14.9 in
Seat height	18.8 / 19.6 / 20.4 / 21.2
Seat depth	13.7 – 21.6 with 0.78 in increments
Footrest	height and angle adjustable
Total weight	46.3 lbs
Batteries	24 Volt, 2.9 Ah Li-Ion

