

# LEVO Summit MA

LEVO - simply straight Standing!  
Innovative Swiss design  
The most active manual – electric  
standing wheelchair  
LEVO Summit EL

LEVO Summit MA



LEVO: Simply straight Standing

LEVO<sup>®</sup>.dt   
THE EXPERTS IN STANDING

# Lightweight standing performance

## Inspired Swiss design

### The first truly active standing wheelchair



The "MA" standing for 100% manual, the LEVO Summit MA is the manual version of the LEVO Summit Series.

The **LEVO Summit MA** belongs worldwide to the most active and lightest Stand-up Wheelchairs with mechanized standing assistance.

The proven standing system allows you to go from seating to standing, or any position in-between. Thanks to its small dimensions, **LEVO Summit MA** is the perfect companion for indoors and outdoors mobility.

#### Highlights:

- Weight from 21 kg
- Weight capacity up to 125 kg
- Automatic backrest compensation "ABC" for ergonomically-correct standing
- Easily configurable and convertible to individual needs
- Easy transfer
- Easy transport

#### Technical Data

Model	S - L
Seat width (SW)	35 / 38 / 40 / 42 / 45 / 48 / 50
Total width (W)	Seat width + 20 cm
Total length (L)	Seat depth + 43 cm
Total height (H) min.	Seat height + 30 cm
Seat height	48 / 50 / 52 / 54
Seat depth	35 – 55 in 2 cm increments
Backrest height	32 – 38 cm
Footrest	height and angle adjustable
Armrest height	22 cm
Total weight minimum	21 kg
Weight capacity	125 kg
Transport volume	Variable W x L x 70...73 cm

#### Standard configuration

M/M Wheelchair with manual drive and manual standing system, flip-up footrests, adjustable knee support

#### Options / Accessories

Seat cushions, push handles, special wheels and many more on request.

#### Colors

Frame in white or black and Details in several colors as standard. Colors according to your wish are on option.

Distributor:

Manufactured by:

LEVO AG | Anglikerstrasse 20 | Tel +41 (0)56 618 44 11 | office@levo.ch | UID-Nr. CHE-107.867.006 MWST  
Switzerland | CH-5610 Wohlen | Fax +41 (0)56 618 44 10 | www.levo.ch | EU SRN CH-MF-000008751

**LEVO**.dt  
THE EXPERTS IN STANDING