

LEVO Dealer:

Date:

Dealer reference / Client name:

LEVO AG

Anglikerstrasse 20

CH-5610 Wohlen

Tel.: +41 56 618 44 11

Fax: +41 56 618 44 10

Office@levo.ch

Order number:


Standard configuration

60.001.7201

- Power Base: width 63 cm; length 105 cm (incl. footplate)
- E/E Wheelchair with 4 active wheels and "Low shearing system"
- Seat support plate
- Seat height 48 cm
- V-Trak Axxis low back, 40 cm/16" high
- Quick release for reclining back
- Flip-up armrests
- VR2 control system from PG Drives Technology
- 48 Ah Batteries and battery charger
- Standard colors: red, yellow, blue, black or white
- Flip-up footrest side- and backwards with non slip surface
- Chest strap and knee support pro
- 8" single rear wheel castors
- Vehicle tie-down locations



Anthropometric Data



Please measure the anthropometric data as shown. If possible, the user should be in a sitting position and wearing shoes.

1 Hip width

Please check the correct box according to the measured hip width in a sitting position:

32 cm ☐ 52 cm

Seat width:

32 cm	☐	Option
36 cm	☐	Option
40 cm	☐	Option
44 cm	☐	Option
48 cm	☐	Option
52 cm	☐	Option

2 Upper leg length

Please fill in the measured upper leg length:

38 cm ☐ 66 cm

3 Lower leg length

Please fill in the measured lower leg length:

24 cm ☐ 44 cm

If the lower leg length is between 45 and 50 cm, C3 is available with a 6 cm higher seat height. The seat height is not adjustable. Please check the box accordingly.

☐ Option

4 Elbow height

Standard armrest height is 22 cm

17 cm ☐ 26 cm

User weight

The weight of the user determines the force of the gas springs. The maximum user weight is 140 kilograms. Please fill in the box.

0 - 60 kg **1200 N**
 60 - 100 kg **1700 N** ☐
 100 - 140 kg **2100 N** kg

Controls



1 VR2 Control system:

- Standard control system
- 90 A Power module
- Max. 2 actuators

- Controls lights
- Programmable

☐ Standard

2 R-Net control system "Standard":

- Color LCD Joystickmodule
- 120 A Power module
- Max. 6 actuators, multiple axes

- Controls lights
- Fully programmable, multiple profiles
- Expandable

☐ 60.072.0010

3 R-Net control system "XL", as R-net control system "Standard" with following extras :



- 3.5" VGA color display
- Paddle switches for on/off, profile/mode and speed

- Metal case
- Loudspeaker instead of a buzzer for a louder horn


☐ 60.072.0020

3 R-Net control system "XL", with integrated bluetooth

☐ 60.072.0030

		Joystick position		Right		Standard	●
				Left		Option	
		① Swing-away, <i>offers more space for transfers and allows to drive closer to tables</i>				60.005.7210	
		② Joystick protection guard, <i>shock protection for joystick, protects the users hands</i>				60.012.6860	
Attendant Controls							
		VR2 Attendant control				32.072.0140	
		R-Net Attendant control				32.068.1850	
		For special controls, see our special controls order form.					
Speed							
		6 km/h				Option	●
		8 km/h				Option	
		10 km/h				Standard	
Stability*							
		Gyro module, <i>only available for R-net control system</i>				32.072.0400	
		<i>- Perfect control on cambers - Straight line driving over obstacles - Eliminates spin-out and fish-tailing</i>					
		<i>* This option is only available in Europe</i>					
Batteries							
		48 Ah				Standard	●
		73 Ah				60.006.7210	
		Without batteries, <i>be aware net price is lower after discount</i>				60.106.6810	
		Without charger, <i>be aware net price is lower after discount</i>				60.111.6710	
Frame color							
		Red	RAL 3020			Standard	●
		Yellow	RAL 1018			Option	
		Blue	RAL 5002			Option	
		Black	RAL 9017			Option	
		White	RAL 9016			Option	
		Blue metallic	RAL 5005-M			60.008.0020	
		Red metallic	RAL 3002-M			60.008.0020	
		Silver metallic	RAL 9007			60.008.0020	
		Orange reflex	RAL 2011-R			60.008.0030	
		Lemon reflex	RAL 6018-R			60.008.0030	
		Silver reflex	RAL 9007-R			60.008.0030	
		Other RAL color	RAL: _____			60.008.0010	
Tube color							
		Standard grey metallic				Standard	●
		Matt black				Option	
Wheels							
		① Block profile drive wheels				Standard	●
		② V-profile drive wheels				Option	
		Puncture proof tires on drive wheels				Standard	●
		Air tires on drive wheels				Option	
		Black tires				On request	
		③ 7" twin rear castors				60.002.0140	
Lights							
			Lights and indicators for outdoor use			60.012.7220	
Seating							
		① Tilt in space 0 - 30°				60.012.7210	
		② Power reclining backrest, <i>only in combination with R-Net, comes with plastic cover</i>				60.009.0150	
		③ Power Elevating Legrests, <i>only in combination with R-Net. Knee support "Integral" fitted.</i> <i>- Length compensated 155mm</i> <i>- Legrests can be controlled separately, but also together</i>				60.012.7230	

Backrest, *The V-Trak Axis low back is delivered with the chair as a standard feature*



- ① V-Trak Axis low backrest, height 40 cm/16"
 ② V-Trak Axis high backrest, height 61 cm/24"
 V-Trak Axis reduced height backrest, height 31 cm/12"

	Standard	•
	60.009.0080	
	60.009.0070	

- ③ V-Trak Axis Deep contour

Size ▶	M	L	XL	XXL		
Type ▼	inside 26 - 37 cm	inside 30 - 40 cm	inside 33 - 46 cm	inside 35 - 52 cm		
a High					60.009.0110	
b Reduced height					60.009.0130	

- ④ Third party backrest conversion bracket, *please choose the correct size*

32 - 36	32.072.0351
40 - 44	32.072.0352
44 - 52	32.072.0353

Without backrest shell, *be aware net price is lower after discount*

	60.109.7210
--	-------------

Seat cushion, *with memory foam depending on the seat depth and seat width incl. cover*



Seat width ▶	32	36	40	44	48	52
Seat depth ▼						
35, 37, 39 cm	30.103.5531	30.103.5532	30.903.5535	30.903.5531	30.903.5538	on request
41, 43, 45 cm	30.903.5532	30.103.5533	30.103.5534	30.103.5539	30.903.5536	30.903.5537
47, 49, 51 cm	30.903.5533	30.103.5535	30.103.5536	30.103.5541	30.103.5542	30.103.5549
53, 55, 57 cm	on request	30.903.5534	30.103.5537	30.103.5543	30.103.5544	30.103.5545
59, 61, 63 cm	on request	on request	30.103.5538	30.103.5546	30.103.5547	30.103.5548

Seat cushion cover

Seat width ▶	32	36	40	44	48	52
Seat depth ▼						
35, 37, 39 cm	30.103.6211	30.103.6212	30.903.6215	30.103.6220	on request	on request
41, 43, 45 cm	30.903.6212	30.103.6213	30.103.6214	30.103.6219	30.103.6235	30.103.6230
47, 49, 51 cm	30.903.6213	30.103.6215	30.103.6216	30.103.6221	30.103.6222	30.103.6229
53, 55, 57 cm	on request	30.903.6214	30.103.6217	30.103.6223	30.103.6224	30.103.6225
59, 61, 63 cm	on request	on request	30.103.6218	30.103.6226	30.103.6227	30.103.6228

Seat cushion cover urin proof for all seat depths and widths

	30.903.6230
--	-------------

Headrest



- ① V-Trak headrest standard, *only in combination with V-Trak Axis backrest*
 ② V-Trak headrest winged, *only in combination with V-Trak Axis backrest*
 Track extension for V-Trak head support

	02.067.116
	02.067.117
	02.009.0028

Knee support



- ① Knee support Pro

Size S (10 - 12 cm)	Standard	•
Size M (14 - 16 cm)	Standard	•

- ② Knee support Integral

Size S (10 - 12 cm)	60.011.0040
Size M (14 - 16 cm)	60.011.0040

Extension for Integral +100 mm

	Option
--	--------

Accessories



- ① Skirt guard

Right	32.072.0021
Left	32.072.0022

- ② Upper leg support

Right	32.072.0031
Left	32.072.0032

Holder to skirt guard and upper leg support

32	36	40	44	48	52
2x 32.072.0151	2x 32.072.0152	2x 32.072.0153	2x 32.072.0154	2x 32.072.0155	2x 32.072.0156

Accessories

③ Chest role

32	36	40	44	48	52
320.720.041	32.072.0042	32.072.0043	32.072.0044	32.072.0045	32.072.0046

④ Foot guide

Small	Right		32.069.95
Small	Left		32.069.96
Large	Right		32.070.0040
Large	Left		32.070.0041

⑤ Push handle bar

	32.072.0090
--	-------------

⑥ Rearview mirror

Right		32.072.0050
Left		32.072.0050

⑦ Swing-away table tray

32 - 40		32.072.0011
44 - 52		32.072.0012

Mounted on the opposite side of the joystick position

Right	
Left	

Quick release kit for seat depth and footrest position adjustment

	60.004.7200
--	-------------

⑧ Joystick top

T-bar		20.067.22
Stick		20.067.23
Softball		20.067.24
Mushroom		20.067.25
Chin		20.067.26
Ball		20.067.27

Remarks