

LEVO Dealer:

Date:

Dealer reference / Client name:

LEVO AG  
 Anglikerstrasse 20  
 CH-5610 Wohlen  
 Tel.: +41 56 618 44 11  
 Fax: +41 56 618 44 10  
 Office@levo.ch

Order number:


Standard configuration

60.001.6801



- Power Base: width 66 cm, length 110 cm (incl. footplate)
- E/E Wheelchair: powered driving, sitting and standing
- Seat plate and backrest upholstery
- Flip-up armrests
- Height adjustable armrests
- Backrest length adjustment and armrests height adjustment
- Manual seat height adjustment
- R-Net control system with color display
- 48 Ah batteries and charger
- Standard colors: white, red, blue or black
- Flip-up footrest
- Chest strap and knee support with wedge

### Anthropometric Data



#### 1 Hip width

Please check the correct box according to the measured hip width in a sitting position:

33 cm  50 cm

##### Seat width:

38 cm	<input type="checkbox"/>	Option
44 cm	<input type="checkbox"/>	Option
50 cm	<input type="checkbox"/>	Option

#### 2 Upper leg length

Please fill in the measured upper leg length, seat depth is 3 cm shorter.

35 cm  59 cm

##### Seat depth:

S-M Frame	36.0 cm	<input type="checkbox"/>	Option
	38.2 cm	<input type="checkbox"/>	Option
	40.4 cm	<input type="checkbox"/>	Option
	42.6 cm	<input type="checkbox"/>	Option
	44.8 cm	<input type="checkbox"/>	Option
L-XL Frame	47.0 cm	<input type="checkbox"/>	Option
	49.2 cm	<input type="checkbox"/>	Option
	51.4 cm	<input type="checkbox"/>	Option
	53.6 cm	<input type="checkbox"/>	Option
	55.8 cm	<input type="checkbox"/>	Option
	58.0 cm	<input type="checkbox"/>	Option

#### 3 Lower leg length

Please fill in the measured lower leg length:

37 cm  64 cm

The seat height is dependent on the lower leg length. With a longer lower leg length the seat is higher.

Lower leg length in cm	37,0-44,2	46,0-47,8	49,6-51,4	53,2-55,0	56,8-58,6	60,4-62,2	64,0
Seat height in cm	56	58.5	61	63.5	66	68.5	71

#### 4 Elbow height

Standard armrest height is 20 cm

17 cm  26 cm

#### Backrest height, for standard backrest

Please measure the distance between seat plate and lowest shoulder blade point:

cm

#### User weight

The weight of the user determines the force of the support wheel gas springs.  
 The maximum user weight is 120 kilograms. Please fill in the box.

kg

Please measure the anthropometric data as shown. If possible, the user should be in a sitting position and wearing shoes.

### Controls



①	R-Net control system "Standard": - Color LCD Joystickmodule - 120 A Power module - Max. 6 actuators, multiple axes - Controls lights - Fully programmable, multiple profiles - Expandable		Standard	●
②	R-Net control system "XL", as R-net control system "Standard" with following extras : - 3.5" VGA color display - Paddle switches for on/off, profile/mode and speed - Metal case - Loudspeaker instead of a buzzer for a louder horn		60.075.0010	
②	R-Net control system "XL" with integrated bluetooth		60.075.0030	

### Joystick position



			Right Left	Standard Option	●
①	Swing-away, offers more space for transfers and allows to drive closer to tables			60.005.6810	
②	Joystick protection guard, shock protection for joystick, protects the users hands			60.012.6860	

### Attendant Controls



	R-Net Attendant control			32.068.1850	
	Connector Block for R-Net, needed for special controls			32.068.1880	
	<i>For special controls, see our special controls order form.</i>				

### Speed

	6 km/h			Option	
	8 km/h			Option	
	10 km/h			Standard	●

### Stability\*



	Gyro module			32.068.1000	
	- Perfect control on cambers - Straight line driving over obstacles - Eliminates spin-out and fish-tailing				
	* This option is only available in Europe				

### Batteries

	48 Ah			Standard	●
	Without batteries, be aware net price is lower after discount				
	Without charger, be aware net price is lower after discount				
				60.106.6810	
				60.111.6710	

### Frame color

White	RAL 9016			Standard	●
Red	RAL 3020			Option	
Blue	RAL 5002			Option	
Black	RAL 9017			Option	
Blue metallic	RAL 5005 - M			60.008.0020	
Red metallic	RAL 3002 - M			60.008.0020	
Silver metallic	RAL 9007			60.008.0020	
Orange reflex	RAL 2011 - R			60.008.0030	
Lemon reflex	RAL 6018 - R			60.008.0030	
Silver reflex	RAL 9007 - R			60.008.0030	
Other RAL color	RAL: _____			60.008.0010	

### Wheels



①	Block profile drive wheels			Standard	●
②	V-profile drive wheels			Option	
	Puncture proof tires on drive wheels			Standard	●
	Air tires on drive wheels			Option	
	Black tires			Option	
③	Curb Climber			60.002.0100	

### Lights



Lights and indicators for outdoor use

60.012.6810

### Seating

### Electrical adjustments



Powered seat lift 150 mm

60.012.6820

① Powered leg rest, *makes an angle adjustment of 90-180° possible, independent from the backrest*

60.012.6870

② Combination of powered reclining back and powered leg rest  
*Enables a complete lying position and the function lying to standing is integrated. Enables a lying transfer.*

60.012.6830

*Recommendation: order with additional V-Trak arms and head rest.*

### Backrests

Standard backrest upholstery

Standard

V-Trak Axxis reduced height backrest, height 31 cm

60.009.0070

V-Trak Axxis low backrest, height 40 cm

60.009.0050

V-Trak Axxis high backrest, height 61 cm

60.009.0060

③ V-Trak Axxis Deep contour

Size	M	L	XL	XXL
Type	inside 26 - 37 cm	inside 30 - 40 cm	inside 33 - 46 cm	inside 35 - 52 cm
a High				
b Reduced height				

60.009.0120

60.009.0140

Without backrest, *b e aware net price is lower after discount*

60.109.6810

### Backrest options

① Backrest extension 15 cm, *for standard backrest*

38

32.062.3815

44

32.062.4415

50

32.062.5015

② Backrest extension with push handle bar and upholstery, *for standard backrest*

Seat width	38	44	50
Seat depth			
15 cm	32.068.0761	32.068.0771	32.068.0781
25 cm	32.068.0762	32.068.0772	32.068.0782
35 cm	32.068.0763	32.068.0773	32.068.0783

Additional V-Trak arms, *recommended when selecting reclining back*

60.009.6810

### Seat cushion with memory foam depending on the seat depth and seat width



Seat width	38	44	50
Seat depth			
S	30.102.9273	30.102.9274	30.102.9275
M	30.102.9276	30.102.9277	30.102.9278
L	30.102.9281	30.102.9282	30.102.9283
XL	30.102.9284	30.102.9285	30.102.9286

### Additional seat cushion cover

Seat width	Standard			Urin Proof		
	38	44	50	38	44	50
Seat depth						
S	30.102.9461	30.102.9462	30.102.9463	30.102.9481	30.102.9482	30.102.9483
M	30.102.9464	30.102.9465	30.102.9466	30.102.9484	30.102.9485	30.102.9486
L	30.102.9467	30.102.9468	30.102.9469	30.102.9487	30.102.9488	30.102.9489
XL	30.102.9471	30.102.9472	30.102.9473	30.102.9474	30.102.9475	30.102.9476

### Headrest



① V-Trak headrest standard, *only in combination with V-Trak and Axxis backrest*

02.067.116

② V-Trak headrest winged, *only in combination with V-Trak and Axxis backrest*  
Track extension for V-Trak head support, *for V-Trak 60 and 80 cm*

02.067.117

02.009.0028

### Lateral supports and hardware



V-Trak lateral support, pad and cover

Medium	Right		32.068.0030
Medium	Left		32.068.0040

V-Trak lateral support, pad and cover

Large	Right		32.068.0010
Large	Left		32.068.0020

### Knee support



Standard knee support

	Standard		
--	----------	--	--

① Rectangular knee wedge

50 x 100		Option	
60 x 100		Option	
80 x 100		Option	
50 x 150		Option	
80 x 150		Option	

② Conical knee wedge, *gives more support, in particular with spasticity*

50-100 x 100		60.011.0020	
90-140 x 100		60.011.0020	
130-180 x 100		60.011.0020	

③ Knee support "Pro"

Foam size S (10 - 12 cm)		60.011.0030	
Foam size M (14 - 16 cm)		60.011.0030	

### Foot plate

Foot plate depth

27 cm		Standard	
20 cm		60.003.0090	

Foot plate width, *30 cm only possible with lower leg lengths between 23 and 27 cm*

22 cm		Standard	
30 cm		60.003.0110	

Foot plate angle, *foot plate is angle adjustable in all three position s*

80°		Option	
90°		Standard	
100°		Option	

① Foot plate wedge, *better positioning with spasticity or foot deformation*

50-100 x 80		30.101.8171	
80-130 x 80		30.101.8172	

② Foot protection frame, *protects the users feet*

		60.003.0100	
--	--	-------------	--

Non-slip footrest surface


		02.001.0702	
--	--	-------------	--

### Accessories



① Skirt guard

S - M		32.068.0830	
L - XL		32.068.0840	



② Upper leg support

Right		32.068.1751	
Left		32.068.1752	



③ Chest role

<b>38</b>	<b>44</b>	<b>50</b>
32.068.1761	32.068.1762	32.068.1763

④ Foot guide

Small	Right		32.069.95
Small	Left		32.069.96
Large	Right		32.070.0040
Large	Left		32.070.0041

⑤ Push handle bar, *can be used with Standard and V-Trak backrest*

<b>38</b>	<b>44</b>	<b>50</b>
32.068.0791	32.068.0792	32.068.0793

⑥ Rearview mirror

Right		32.068.0120	
Left		32.068.0120	

⑦ Swing-away table tray, *can only be mounted on the opposite side of the joystick position*

<b>38</b>	<b>44</b>	<b>50</b>
32.068.1061	32.068.1062	32.068.1062

Right		
Left		

**Accessories**



⑧ Fixed table tray

<b>38</b>	<b>44</b>	<b>50</b>
32.062.0038	32.062.0044	32.062.0050



⑨ Support for QLK docking system, *without assembly cost*  
 Assembly cost for EZ lock brackets

Vehicle tie-down kit	32.068.1612
	32.068.2050
	33.000.0000

**Accessories**



⑩ Joystick top

<b>T-bar</b>	20.067.22
<b>Stick</b>	20.067.23
<b>Softball</b>	20.067.24
<b>Mushroom</b>	20.067.25
<b>Chin</b>	20.067.26
<b>Ball</b>	20.067.27

Remarks